



REVISION SCHEDULE		
REVISION #	DESCRIPTION	PLANNER/DATE




## STERILE PROCESSING DEPARTMENT FLOOR PLAN

**\*THIS SAMPLE LAYOUT IS INTENDED TO ENCOURAGE DISCUSSIONS FOR AN INTEGRATED DESIGN PROCESS. FOR A COMPREHENSIVE CAPACITY ASSESSMENT AND OPTIMIZED DESIGN FOR YOUR PROCESSING DEPARTMENT, PLEASE CONTACT YOUR LOCAL STERIS REPRESENTATIVE.**

Processing Capacity up to 300 Trays Hours Operation - 24 Peak Processing Hours - 10		
FLOOR PLAN Design Template Sterile Processing Department PROPOSED EQUIPMENT LAYOUT		
 SYSTEMS PROJECT PLANNING STERIS Corporation 5900 Healey Road Mentor, OH 44860 USA 800-548-4837	Dr: t11	H-
	Ck:	D-
Date: 23 October 2013	Scale: 3/16" = 1'-0"	PAGE 1 OF 3

REVISION SCHEDULE		
REVISION #	DESCRIPTION	PLANNER/DATE



FLOW LEGEND	
CASE CARTS	→
INSTRUMENTS	→
HAND WASH ITEMS	→
LOW TEMP INSTRUMENTS	→
UTENSIL CARTS	→
STAFF	→

## STERILE PROCESSING DEPARTMENTAL FLOW

**\*THIS SAMPLE LAYOUT IS INTENDED TO ENCOURAGE DISCUSSIONS FOR AN INTEGRATED DESIGN PROCESS. FOR A COMPREHENSIVE CAPACITY ASSESSMENT AND OPTIMIZED DESIGN FOR YOUR PROCESSING DEPARTMENT, PLEASE CONTACT YOUR LOCAL STERIS REPRESENTATIVE.**

Processing Capacity up to 300 Trays Hours Operation - 24 Peak Processing Hours - 10		
FLOOR PLAN Design Template Sterile Processing Department PROPOSED EQUIPMENT LAYOUT		
 <small>STERIS CORPORATION          5900 Healey Road          Mentor, OH 44860 USA          800-548-4837</small>	Dr: t11 Ck:	H- D-
	Date: 23 October 2013 Scale: 3/16" = 1'-0"	PAGE: 1 OF 3

REVISION SCHEDULE		
REVISION #	DESCRIPTION	PLANNER/DATE



STERILE PROCESSING OVERALL AREA: 6000-6700 ft <sup>2</sup>	(100%)
DECONTAMINATION AREA: 1050-1150 ft <sup>2</sup>	(~20%)
ASSEMBLY AREA: 1950-2050 ft <sup>2</sup>	(~30%)
STERILIZATION AREA: 700-800 ft <sup>2</sup>	(~10%)
PROCESSED STORAGE & SUPPLY AREA: 1000-1200 ft <sup>2</sup>	(~20%)
ADMINISTRATION & SUPPORT AREA: 1300-1500 ft <sup>2</sup>	(~20%)

## STERILE PROCESSING DEPARTMENT AREA DISTRIBUTION

Processing Capacity up to 300 Trays Hours Operation - 24 Peak Processing Hours - 10		
FLOOR PLAN Design Template Sterile Processing Department PROPOSED EQUIPMENT LAYOUT		
 <small>STERIS CORPORATION 3900 Healey Road Meriden, CT 06450 USA 800-548-4837</small>	Dr: t11	H-
	Ck:	D-
Date: 23 October 2013	Scale: 3/16" = 1'-0"	PAGE 1 OF 3

**\*THIS SAMPLE LAYOUT IS INTENDED TO ENCOURAGE DISCUSSIONS FOR AN INTEGRATED DESIGN PROCESS. FOR A COMPREHENSIVE CAPACITY ASSESSMENT AND OPTIMIZED DESIGN FOR YOUR PROCESSING DEPARTMENT, PLEASE CONTACT YOUR LOCAL STERIS REPRESENTATIVE.**